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- 1 ON/OFF Button
- 2 Display
- 3 Slot-in Card
- 4 Cuff Socket
- 5 Mains Adapter Socket
- 6 Cuff Compartment
- 7 Battery Compartment
- 8 Cuff
- 9 Cuff Connector
- AT M-Button (Memory)
- ⑪ Time Button

Display

- ⑫ Heart Arrhythmia Indicator
- AM Pulse Rate
- AN Battery Display
- AO Stored Value
- AP Systolic Value
- AQ Diastolic Value
- AR Pulse
- AS Date/Time
- BT Radio clock symbol

Dear Customer,

Your new aponorm® by microlife blood pressure monitor is a reliable medical instrument for taking measurements on the upper arm. It is simple to use, accurate and comes highly recommended for blood pressure monitoring in your home. This instrument was developed in collaboration with physicians and clinical tests prove its measurement accuracy to be very high.

Please read through these instructions carefully so that you understand all functions and safety information. We want you to be happy with your aponorm® by microlife product. If you have any questions, problems or want to order spare parts, please contact aponorm® by microlife-Customer Service. Your pharmacy will be able to give you the address of the aponorm® by microlife dealer in your country. Alternatively, visit the Internet at www.aponorm.de where you will find a wealth of invaluable information on our products.

Stay healthy – aponorm® by microlife!

Table of Contents

1. Important Facts about Blood Pressure and Self-Measurement

- How do I evaluate my blood pressure?

2. Using the Instrument for the first Time

- Activate the fitted batteries
- Setting the date and time
- Switching beep On/Off
- Select the correct cuff

3. Taking a Blood Pressure Measurement using this Instrument

4. Appearance of the Heart Arrhythmia Indicator for early Detection

5. Data Memory

- Viewing the stored values
- Memory full
- Delete a measurement
- Clear all values

6. Replacing the Slot-in Card

7. Battery Indicator and Battery change

- Batteries almost flat
- Batteries flat – replacement
- Which batteries and which procedure?
- Using rechargeable batteries

8. Using a Mains Adapter

9. Error Messages

10. Safety, Care, Accuracy Test and Disposal

- Safety and protection
- Instrument care
- Cleaning the cuff
- Accuracy test
- Disposal

11. Guarantee

12. Technical Specifications

Guarantee Card (see Back Cover)

1. Important Facts about Blood Pressure and Self-Measurement

- **Blood pressure** is the pressure of the blood flowing in the arteries generated by the pumping of the heart. Two values, the **systolic** (upper) value and the **diastolic** (lower) value, are always measured.
- The instrument also indicates the **pulse rate** (the number of times the heart beats in a minute).
- **Permanently high blood pressure values can damage your health and must be treated by your doctor!**
- Always discuss your values with your doctor and tell him/her if you have noticed anything unusual or feel unsure. **Never rely on single blood pressure readings.**
- Enter your readings in the enclosed **blood pressure diary**. This will give your doctor a quick overview.
- There are many causes of excessively **high blood pressure values**. Your doctor will explain them in more detail and offer treatment where appropriate. Besides medication, relaxation techniques, weight loss and exercise can also lower your blood pressure.
- **Under no circumstances should you alter the dosages of any drugs prescribed by your doctor!**
- Depending on physical exertion and condition, blood pressure is subject to wide fluctuations as the day progresses. **You should therefore take your measurements in the same quiet conditions and when you feel relaxed!** Take at least two measurements per day, one in the morning and one in the evening.
- It is quite normal for two measurements taken in quick succession to produce significantly **different results**.
- **Deviations** between measurements taken by your doctor or in the pharmacy and those taken at home are quite normal, as these situations are completely different.
- **Several measurements** provide a much clearer picture than just one single measurement.
- **Leave a small break** of at least 15 seconds between two measurements.
- If you are **pregnant**, you should monitor your blood pressure very closely as it can change drastically during this time!

- If you suffer from an **irregular heartbeat** (arrhythmia, see «Section 4.»), measurements taken with this instrument should only be evaluated after consultation with your doctor.
- The pulse display is not suitable for checking the frequency of heart pacemakers!**

How do I evaluate my blood pressure?

Table for classifying blood pressure values in adults in accordance with the World Health Organisation (WHO) in 2003. Data in mmHg.

Range	Systolic	Diastolic	Recommendation
blood pressure too low	↓ 100	↓ 60	Consult your doctor
1. blood pressure optimum	100 - 120	60 - 80	Self-check
2. blood pressure normal	120 - 130	80 - 85	Self-check
3. blood pressure slightly high	130 - 140	85 - 90	Consult your doctor
4. blood pressure too high	140 - 160	90 - 100	Seek medical advice
5. blood pressure far too high	160 - 180	100 - 110	Seek medical advice
6. blood pressure dangerously high	180 ↑	110 ↑	Urgently seek medical advice!

The higher value is the one that determines the evaluation.

Example: a readout value between **150/85** or **120/98** mmHg indicates «blood pressure too high».

The slot-in card **3** on the front of the instrument shows ranges 1-6 in the Table.

2. Using the Instrument for the First Time

Activate the fitted batteries

Pull out the protective strip projecting from the battery compartment **7**.

Setting the date and time

Your aponorm® Basis Plus upper arm blood pressure monitor is equipped with a radio clock, which is controlled by the DCF77 signal. This time and date are set automatically. If no DCF77 signal (this is the term for the time signal transmitter) is available, the time is not adjusted automatically.

The display comes to receiving mode and device sets the time automatically, if you drag the zip to activate the batteries (or insert batteries).

☞ If you do not wish to set the time or you want to take a measurement immediately, exit the time setting mode by pressing the On/Off button.

Once the time signal is received, the radio clock icon  BT flashes on the display every second until the reception of the DCF77 signal is complete.

While receiving the DCF77 signal, the radio clock icon  BT flashes first with 1, then 2 and then 3 bars above it. When the icon appears with 3 bars, the time has been set. First the date and then the time are displayed at the top right of the display. Then the time is displayed permanently together with the radio clock icon  BT.

☞ The automatic setting of date and time takes about 2-4 minutes. If no signal is received within 10 minutes, „- - -“ appears on the display.

Now you need to improve the reception by choosing another location or set the clock manually (for more information see „manual setting“). After a change of location, batteries have to be re-inserted or to reset radio clock to „ON“ to restart the auto setting again.

Manual setting (if automatic setting does not work or is not desired):

The device must be in standby mode. If you do not press any buttons for 1 minute during the time setting, the time setting mode switches off and „- : - -“ appears in the clock display.

1. Press down Time button several seconds until „bl OFF“ appears on the display.
2. Press Time button again, the radio clock  indicator and „ON“ appears in the display.
3. Press the M button to switch on („ON“) or switch off the radio clock („OFF“).
4. Press ON/OFF button to confirm and return to standby mode, or if the radio clock is switched off, you can press Time button to confirm (continue with step 5).
5. The display shows a flashing year.
 - a. You can press the M button to set the year.
 - b. Pressing the clock button flashes the month. You can press the M button to set the month.
 - c. Pressing the clock button flashes the day. You can press the M button to set the day.
 - d. Pressing the clock button flashes the hour. You can press the M key to set the hour.
 - e. Pressing the clock button flashes the minute. You can press the M button to set the minute.
 - f. Please press once again the clock button to complete the action. Then the date and time you set are shown twice each in succession before the device returns to standby mode. In standby mode, the time is displayed permanently.

 If you press the On/Off button during the manual time setting, you will exit the time setting mode. The so far entered date and time values are kept. „- : - -“ appears in the display if you did not set the time previously.

Switching beep On/Off

There is no beeps sound when you take a measurement using the device. You can switch on the beeps by making the following setting.

To switch between „beep off“ and „beep on“, press and hold the Time-button until bl and OFF appear flashing in the display.

Press the M button to switch between OFF and ON.

Press the Time button, device shows the next setting: „Switching radio clock on/off“. You can then press On/Off button to confirm the beep on/off selection or go with the radio clock on/off settings. If you select beep „OFF“, the beep function remains switched off when you take measurements.

Select the correct cuff

aponorm® by microlife offers 3 different cuff sizes: S, M and L. Select the cuff size to match the circumference of your upper arm (measured by close fitting in the centre of the upper arm). M-L is the correct size for most people.

Cuff size	for circumference of upper arm
M-L	22 - 42 cm (8.75 - 16.5 inches)
S	17 - 22 cm (6.75 - 8.75 inches)
M	22 - 32 cm (8.75 - 12.5 inches)
L	32 - 42 cm (12.5 - 16.5 inches)

 Use only aponorm® by microlife cuffs!

- ▶ Contact aponorm® by microlife Service, if the enclosed cuff ⑧ does not fit.
- ▶ Connect the cuff to the instrument by inserting the cuff connector ♀ into the cuff socket 4 as far as it will go.

3. Taking a Blood Pressure Measurement using this Instrument

Checklist for taking a reliable measurement

1. Avoid activity, eating or smoking immediately before the measurement.
2. Sit down for at least 5 minutes before the measurement - and relax.
3. Always measure on the same arm (normally left).
4. Remove close-fitting garments from the upper arm. To avoid constriction, shirt sleeves should not be rolled up - they do not interfere with the cuff if they are laid flat.
5. Always ensure the cuff is positioned correctly, as shown in the pictures appearing at the beginning of this booklet.
 - Fit the cuff closely, but not too tight.
 - Make sure that the cuff is 3 cm (1 inch) above your elbow with the tube on the inside of your arm.

- Support your arm so it is relaxed.
 - Ensure that the cuff is at the same height as your heart.
6. Press the On/Off button 1 to start the measurement.
7. The cuff will now pump up automatically. Relax, do not move and do not tense your arm muscles until the measurement result is displayed. Breathe normally and do not talk.
8. When the correct pressure is reached, the pumping stops and the pressure falls gradually. If the required pressure was not reached, the instrument will automatically pump some more air to the cuff.
9. During the measurement, the heart icon 13 is flashing.
10. The result, comprising the systolic 16 and the diastolic 17 blood pressure and the pulse 18 is displayed and longer beep is heard. Note also the explanations on further displays in this booklet.
11. When the measurement has finished, remove the cuff and pack it into the instrument as shown in Fig. II.
12. Enter the result in the enclosed blood pressure pass and switch off the instrument. (The monitor does switch off automatically after approx. 1 min.).

 You can stop the measurement at any time by pressing the ON/OFF button (e.g. if you feel uneasy or an unpleasant pressure sensation).

4. Appearance of the Heart Arrhythmia Indicator for early Detection

This symbol 12 indicates that certain pulse irregularities were detected during the measurement. In this case, the result may deviate from your normal blood pressure – repeat the measurement. In most cases, this is no cause for concern. However, if the symbol appears on a regular basis (e.g. several times a week with measurements taken daily) we advise you to tell your doctor.

Please show your doctor the following explanation:

Information for the doctor on frequent appearance of the Arrhythmia indicator

This instrument is an oscillometric blood pressure monitor that also analyses pulse frequency during measurement. The instrument is clinically tested.

The arrhythmia symbol is displayed after the measurement, if pulse irregularities occur during measurement. If the symbol appears more frequently (e.g. several times per week on measurements performed daily) we recommend the patient to seek medical advice.

The instrument does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

5. Data Memory

At the end of a measurement, this instrument automatically stores each result, including date and time.

Viewing the stored values

Press the M-button AT briefly, when the instrument is switched off. The display first shows «M» 15 and then a value, e.g. «M 17». This means that there are 17 values in the memory. The instrument then switches to the last stored result.

Pressing the M-button again displays the previous value. Pressing the M-button repeatedly enables you to toggle between one stored value and another.

Memory full

 When the memory has stored 60 results, the memory is full. From this point onwards, a new measured value is stored by **overwriting the oldest value**.

Delete a measurement

Select the measurement value you want to delete. Press the ON/OFF Button and keep it pressed for approx. 10 seconds until the memory symbol flashes. Confirm your deletion request by pressing the memory button. "CL" appears briefly on the display and the value you selected has been deleted. The display returns to standby mode.

Clear all values

If you are sure that you want to permanently remove all stored values, hold down the M-button (the instrument must have been switched off beforehand) until «CL» appears and then release the button. To permanently clear the memory, press the M-button while «CL» is flashing.

6. Replacing the Slot-in Card

You can replace the slot-in card 3 by pulling it out to the side, as shown in Fig. IV and replacing the paper insert.

It may be helpful to have your doctor note down the medication dosage or an emergency telephone number on the card. Extra cards are supplied with the instrument for this purpose.

7. Battery Indicator and Battery change

Batteries almost flat

When the batteries are approximately $\frac{3}{4}$ used the battery symbol 14 will flash as soon as the instrument is switched on (partly filled battery displayed). Although the instrument will continue to measure reliably, you should obtain replacement batteries.

Batteries flat – replacement

When the batteries are flat, the battery symbol 14 will flash as soon as the instrument is switched on (flat battery displayed). You cannot take any further measurements and must replace the batteries.

1. Open the battery compartment 7 at the back of the instrument by pushing inwards at the two arrows and pulling out the battery compartment cover.
 2. Replace the batteries – ensure correct polarity as shown by the symbols in the compartment.
 3. To set date and time, follow the procedure described in «Section 2.».
- ☞ The memory retains all values although date and time (and possibly also set alarm times) must be reset – the year number therefore flashes automatically after the batteries are replaced.

Which batteries and which procedure?

- ☞ Please use 4 new, long-life 1,5 V, size AA batteries.
- ☞ Do not use batteries beyond their date of expiry.
- ☞ Remove batteries, if the instrument is not going to be used for a prolonged period.

Using rechargeable batteries

You can also operate this instrument using rechargeable batteries.

- ☞ Please use only type «NiMH» reusable batteries!
- ☞ The batteries must be removed and recharged, if the battery symbol (battery flat) appears! They must not remain inside the instrument, as they may become damaged (total discharge as a result of low use of the instrument, even when switched off).
- ☞ Always remove the rechargeable batteries, if you do not intend to use the instrument for a week or more!
- ☞ The batteries can NOT be charged in the blood pressure monitor! Recharge these batteries in an external charger and observe the information regarding charging, care and durability!

8. Using a Mains Adapter

You can operate this instrument using the aponorm® by microlife adapter (DC 6V, 600mA).

- ☞ Use only the aponorm® by microlife mains adapter available as an original accessory appropriate for your supply voltage, e.g. the «aponorm® by microlife 230V adapter».
- ☞ Ensure that neither mains adapter nor cable are not damaged.

 1. Plug the adapter cable into the mains adapter socket 5 in the blood pressure monitor.
 2. Plug the adapter plug into the wall socket.

When the mains adapter is connected, no battery current is consumed.

9. Error Messages

If an error occurs during the measurement, the measurement is interrupted and an error message, e.g. «ERR 3», is displayed.

Error	Description	Potential cause and remedy
«ERR 1»	Signal too weak	The pulse signals on the cuff are too weak. Re-position the cuff and repeat the measurement.*
«ERR 2»	Error signal	During the measurement, error signals were detected by the cuff, caused for instance by movement or muscle tension. Repeat the measurement, keeping your arm still.
«ERR 3»	No pressure in the cuff	An adequate pressure cannot be generated in the cuff. A leak may have occurred. Check that the cuff is correctly connected and is not too loose. Replace the batteries if necessary. Repeat the measurement.
«ERR 5»	Abnormal result	The measuring signals are inaccurate and no result can therefore be displayed. Read through the checklist for performing reliable measurements and then repeat the measurement.*
«HI»	Pulse or cuff pressure too high	The pressure in the cuff is too high (over 300 mmHg) OR the pulse is too high (over 200 beats per minute). Relax for 5 minutes and repeat the measurement.*
«LO»	Pulse too low	The pulse is too low (less than 40 beats per minute). Repeat the measurement.*

* Please consult your doctor, if this or any other problem occurs repeatedly.

 If you think the results are unusual, please read through the information in «Section 1.» carefully.

10. Safety, Care, Accuracy Test and Disposal

Safety and protection

- This instrument may be used only for the purpose described in this booklet. The manufacturer cannot be held liable for damage caused by incorrect application.
- This instrument comprises sensitive components and must be treated with caution. Observe the storage and operating conditions described in the «Technical Specifications» section!
- Protect it from:
 - water and moisture
 - extreme temperatures
 - impact and dropping
 - contamination and dust
 - direct sunlight
 - heat and cold
- The cuffs are sensitive and must be handled with care.
- Only pump up the cuff when fitted.
- Do not use the instrument close to strong electromagnetic fields such as mobile telephones or radio installations.
- Do not use the instrument if you think it is damaged or notice anything unusual.
- Never open the instrument.
- If the instrument is not going to be used for a prolonged period the batteries should be removed.
- Read the further safety instructions in the individual sections of this booklet.

Ensure that children do not use the instrument unsupervised; some parts are small enough to be swallowed.

Instrument care

Clean the instrument only with a soft, dry cloth.

Cleaning the cuff

You can machine wash the **cuff cover** at 30°C (do not iron!).

⚠️ WARNING: Under no circumstances, however, may you wash the inner bladder! Always remove the sensitive bladder from the sleeve before washing and replace it carefully again afterwards.

Accuracy test

We recommend this instrument is tested for accuracy every 2 years or after mechanical impact (e.g. being dropped). Please contact aponorm® by microlife-Service to arrange the test (see foreword).

 Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

11. Guarantee

This instrument is covered by a **5 year guarantee** from the date of purchase. The guarantee is valid only on presentation of the guarantee card completed by the dealer (see back) confirming date of purchase or the till receipt.

- The guarantee covers both instrument and cuff. Batteries are not included.
- Opening or altering the instrument invalidates the guarantee.
- The guarantee does not cover damage caused by improper handling, discharged batteries, accidents or non-compliance with the operating instructions.

Please contact aponorm® by micrlife-Service (see foreword).

12. Technical Specifications

Operating temperature: 10 - 40 °C / 50 - 104 °F

Storage temperature: -20 - +50 °C / -4 - +122 °F

15 - 90 % relative maximum humidity
690 g (including batteries)

Weight:

160 x 120 x 98 mm

Measuring procedure: oscilloscopic, corresponding to Korotkoff method: Phase I systolic, Phase V diastolic

Measurement range: 30 - 280 mmHg – blood pressure
40 - 200 beats per minute – pulse

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GmbH & Co. KG
Microlife AG